

Product Datasheet

Bat IL-8 Recombinant Protein (orb1819699)

Description	The Black flying fox Bat IL-8 yeast-derived recombinant protein is not tagged and is naturally endot...
Target	IL-8
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	AVLSRISTEL RCQCIKTHSA PFHPKFIKEL RVIESGPHCA HTEIIVKLVK GNEVCLNPKE KVVQMVVQAF LKRAEKKDP (79)
Purity	98%
Source	Yeast
MW	9.0 kDa
Entrez	102892534
Expiration Date	6 months from date of receipt.