

Product Datasheet

Ferret CCL1 Recombinant Protein (orb1819710)

Description	The Ferret CCL1 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	CCL1
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	RSMHVSSSTC CFVFVEKRIP QHIIQCYKPT SSSCFHKGLI IKLKGGRETG VSDSAQWVQR YLKTAKSCWV KRT (73)
Purity	98%
Source	Yeast
MW	8.3 kDa
Entrez	101687582
Expiration Date	6 months from date of receipt.