

## Product Datasheet

### Mallard IL-4 Recombinant Protein (orb1819712)

<b>Description</b>	The Duck IL-4 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The p...
<b>Target</b>	IL-4
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	HEAAVLERSI FLKESIRLLG EILSTQVSCD KANVTNVFAG NETGDMELL CKASTVVFES LSCHKPLKGI YLNLLHIVTK STSLKAPCPV AAGNTTSLQE FLRGLHRTLQ RVAKENL (117)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	12.8 kDa
<b>Entrez</b>	<a href="#">101801081</a>
<b>Expiration Date</b>	6 months from date of receipt.