

Product Datasheet

Mallard IL-4 Recombinant Protein (orb1819713)

Description	The Duck IL-4 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The p...
Target	IL-4
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	HEAAVLERSI FLKESIRLLG EILSTQVSCD KANVTNVFAG NETGDMELL CKASTVVFES LSCHKPLKGI YLNLLHIVTK STSLKAPCPV AAGNTTSLQE FLRGLHRTLQ RVAKENL (117)
Purity	98%
Source	Yeast
MW	12.8 kDa
Entrez	101801081
Expiration Date	6 months from date of receipt.