

Product Datasheet

Ferret IL-8 Recombinant Protein (orb1819716)

Description	The Ferret IL-8 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	IL-8
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	EAEAAVLTRV SSELRCQCIK THSTPFHPKY IKELRVIDSG PHCENSEIIV KLVNGKEVCL DPKEDWVQRI VQIFLKKAEK QAA (83)
Purity	98%
Source	Yeast
MW	9.4 kDa
Entrez	101676365
Expiration Date	6 months from date of receipt.