

Product Datasheet

Ferret CCL3 Recombinant Protein (orb1819724)

| | |
|-------------------------|---|
| Description | The Ferret/Mink CCL3 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free... |
| Target | CCL3 |
| Form/Appearance | Lyophilized |
| Storage | -20°C |
| Note | For research use only |
| Protein Sequence | APFGADTPTA CCFSYVSKQI TRKFIVDYFE TSSQCSKPGV IFQTRKGRQL CADPTEDWVQ EYVTDLELNA (70) |
| Purity | 98% |
| Source | Yeast |
| MW | 7.9 kDa |
| Entrez | 101676470 |
| Expiration Date | 6 months from date of receipt. |