

Product Datasheet

Ferret CCL3 Recombinant Protein (orb1819725)

Description	The Ferret/Mink CCL3 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
Target	CCL3
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	APFGADTPTA CCFSYVSKQI TRKFIVDYFE TSSQCSKPGV IFQTRKGRQL CADPTEDWVQ EYVTDLELNA (70)
Purity	98%
Source	Yeast
MW	7.9 kDa
Entrez	101676470
Expiration Date	6 months from date of receipt.