

Product Datasheet

Mouse Midkine Recombinant Protein (orb1819730)

Description	The Mouse Midkine yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. T...
Target	Midkine
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	KKKEKVKKGS ECSEWTWGPC TPSSKDCGMG FREGTCGAQT QRVHCKVPCN WKKEFGADCK YKFESWGACD GSTGTKARQG TLKKARYNAQ CQETIRVTKP CTSKTKSKTK AKKGGKGD (118)
Purity	98%
Source	Yeast
MW	13.1 kDa
Entrez	17242
Expiration Date	6 months from date of receipt.