

## Product Datasheet

### Ferret IL-15 Recombinant Protein (orb1819736)

|                         |  |
|-------------------------|--|
| <b>Description</b>      | The Ferret IL-15 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...  |
| <b>Target</b>           | IL-15  |
| <b>Form/Appearance</b>  | Lyophilized  |
| <b>Storage</b>          | -20°C  |
| <b>Note</b>             | For research use only  |
| <b>Protein Sequence</b> | GIHVFILGSI SAGLPKTEAK WEFVIKDLEK<br>IDKIIQSIHI DTTLYTESDA HPGCKVTAMN<br>CFLLELRVVL LEFSHANLLN ETIHNIIIQA<br>NSNLSSNANI TETGCKECEE LEEKNIKEFF<br>QSFLHIVQIF HNS (133) |
| <b>Purity</b>           | 98%  |
| <b>Source</b>           | Yeast  |
| <b>MW</b>               | 15.0 kDa   |
| <b>Entrez</b>           | <b>101680433</b>   |
| <b>Expiration Date</b>  | 6 months from date of receipt.   |