

Product Datasheet

Ferret IL-15 Recombinant Protein (orb1819737)

Description	The Ferret IL-15 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	IL-15
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	GIHVFILGSI SAGLPKTEAK WEFVIKDLEK IDKIIQSIHI DTTLYTESDA HPGCKVTAMN CFLLELRVVL LEFSHANLLN ETIHNIIIQA NSNLSSNANI TETGCKECEE LEEKNIKEFF QSFLHIVQIF HNS (133)
Purity	98%
Source	Yeast
MW	15.0 kDa
Entrez	101680433
Expiration Date	6 months from date of receipt.