

Product Datasheet

Canine M-CSF Recombinant Protein (orb1819738)

Description	The Canine M-CSF yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	M-CSF
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	EEVSEHCSHM IGNGHLQFLQ QLIDSQMETS CQIAFEFVDQ ELLKDPVCYL KKAFVLMQDI MEDTMRFKDN TPNANVIVKL QELSLRLKSC FTKDYEEQDK ACVRTFYETP LQLLEKIKNV FNETKNLLKK DWNIFSKNCN NSFACSSQG HEKQQEGP (158)
Purity	98%
Source	Yeast
MW	18.5 kDa
Entrez	611795
Expiration Date	6 months from date of receipt.