

Product Datasheet

Bovine CD27 Recombinant Protein (orb1819745)

Description	The Bovine CD27 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	CD27
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	TPGPKSCPEK HYWAQGGWCC QMCEPGTFLV KDCEQHREAA QCNPCTPGVS FMPDHHSRPH CESCRHCNSG LLIRNCTLTA NSECACPEGQ QCRDKDCMEC DGPAQAPGPH PQPSHLPYAE EIAESRTDRH TQTLANSRWL PAPTLSTHWS PQRSLCSANC (160)
Purity	98%
Source	Yeast
MW	17.7 kDa
Entrez	512514
Expiration Date	6 months from date of receipt.