

Product Datasheet

Mouse CXCL4 Recombinant Protein (orb1819757)

Description	The Mouse CXCL4 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	CXCL4
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	VTSAGPEESD GDLSCVCVKT ISSGIHLKHI TSLEVIKAGR HCAVPQLIAT LKNRKCICLD RQAPLYKKVI KKILES (76)
Purity	98%
Source	Yeast
MW	8.2 kDa
Entrez	56744
Expiration Date	6 months from date of receipt.