

Product Datasheet

Human B7-H2 Recombinant Protein (orb1819762)

Description	The Human B7-H2 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	B7-H2
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	DTQEKEVRAM VGSDVELSCA CPEGSRFDLN DVYVYWQTSE SKTVVTYHIP QNSSLENVDS RYRNRALMSP AGMLRGDFSL RLFNVTPQDE QKFHCLVLSQ SLGFQEVLVS EVTLHVAANF SVPVVSAPHS PSQDELTFTC TSINGYPRPN VYWINKTDNS LLDQALQNDT VFLNMRGLYD VVSVLRIART PSVNIGCCIE NVLLQQNLTV GSQTGNDIGE RDKITENPVS TGEKNAAT (238)
Purity	98%
Source	Yeast
MW	26.4 kDa
Entrez	23308
Expiration Date	6 months from date of receipt.