

Product Datasheet

Feline IL-25 Recombinant Protein (orb1819764)

Description	The Feline IL-25 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	IL-25
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	IRIKKDCIHW SSCPSKGQN PTHEWL TENT VLMTTPESAS LIHSLESCRA SEDGPLNSRS ISPWRYELDR DLNRLPQDLY HARCLCPHCV SLQTGSHMDP LGNSELLYHN QTVFYRRPCP GEQGTRDGYC LEQRLYRVSL ACVCVRPRVM A (151)
Purity	98%
Source	Yeast
MW	17.3 kDa
Entrez	101095746
Expiration Date	6 months from date of receipt.