

Product Datasheet

Human IL-1 alpha Recombinant Protein (orb1819765)

Description	The Human IL-1 alpha yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
Target	IL-1 alpha
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SAPFSFLSNV KYNFMRIIKY EFILNDALNQ SIIRANDQYL TAAALHNLDE AVKFDMGAYK SSKDDAKITV ILRISKTKLY VTAQDEDQPV LLKEMPEIPK TITGSETNLL FFWETHGTKN YFTSVAHPNL FIATKQDYWV CLAGGPPSIT DFQILENQ (158)
Purity	98%
Source	Yeast
MW	18.0 kDa
Entrez	3552
Expiration Date	6 months from date of receipt.