

Product Datasheet

Swine IL-16 Recombinant Protein (orb1819769)

Description	The Swine IL-16 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	IL-16
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	ASPDLSSTD SAASVSVASD ISIDSTEATV CTVTLEK TSA GLGFSLEGGK GSLHGDKPLT INRIFRGAAS EQNETIQPGD EILHLAGTAV QGLTRFEAWN IIKALPDGPV TIVIRKSLQ SKGTTGAGDP (130)
Purity	98%
Source	Yeast
MW	13.4 kDa
Entrez	396589
Expiration Date	6 months from date of receipt.