

Product Datasheet

Canine IL-3 Recombinant Protein (orb1819773)

Description	The Canine IL-3 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	IL-3
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	GRPFSTDLPK QYFTMINEIM EMLNKSPSPS EEPLDSNEKE TLEEDTLLRP NLDVFLNASS KFHKNGLLIW>NNLKEFLPLL PTPTRGEPI SIMENNWGDF QRKLKYLEA LDNFLNFKNK P (121)
Purity	98%
Source	Yeast
MW	13.1 kDa
Entrez	481497
Expiration Date	6 months from date of receipt.