

Product Datasheet

Swine IL-5 Recombinant Protein (orb1819776)

Description	The Swine IL-5 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The ...
Target	IL-5
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	VENTMNRLVA ETLTLLSIHR TLLIGDGNLM ISTPVHTNHQ LCIEEVFQGI DTLKNQTARG DAVEKLFQNL SLIKEYIDRQ KKNCGGERWR VTQFLDYLVQ FLGVINTEWT MES (113)
Purity	98%
Source	Yeast
MW	13.0 kDa
Entrez	397409
Expiration Date	6 months from date of receipt.