

Product Datasheet

Swine IL-25 Recombinant Protein (orb1819777)

Description	The Swine IL-25 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	IL-25
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	WSQKGCPHWP NCCPSKGQNP IEEWLKQNAV LMPPLEMASP TPHPESCKAS EDGPLNSRSI APWRYELDRD LNRLPQDLYH ARCLCPHCVS LRTGSHMDPL GNSELLYHNQ TVFYRRPCPG QQGAHDSYCL EPRLYQVSLA CVCVRPRVMA (150)
Purity	98%
Source	Yeast
MW	17.1 kDa
Entrez	100154041
Expiration Date	6 months from date of receipt.