

Product Datasheet

Mouse IL-25 Recombinant Protein (orb1819778)

Description	The Mouse IL-25 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	IL-25
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	VSLRIQEGCS HLPSCCPSKE QEPPEEWLKW SSASVSPPEP LSHTHHAESC RASKDGPLNS RAISPWSYEL DRDLNRVPQD LYHARCLCPH CVSLQTGSHM DPLGNSVPLY HNQTVFYRRP CHGEEGTHRR YCLERRLYRV SLACVCRPR VMA (153)
Purity	98%
Source	Yeast
MW	17.5 kDa
Entrez	140806
Expiration Date	6 months from date of receipt.