

Product Datasheet

Human CXCL3 Recombinant Protein (orb1819784)

Description	The Human CXCL3 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	CXCL3
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	ASVVT EL R RC Q CLQTLQGIHL KNIQSVNVRS PGPHCAQTEV IATLKNGKKA CLNPASPMVQ KIIEKILNKG STN (73)
Purity	98%
Source	Yeast
MW	7.9 kDa
Entrez	2921
Expiration Date	6 months from date of receipt.