

Product Datasheet

Mouse IL-16 Recombinant Protein (orb1819789)

Description	The Mouse IL-16 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	IL-16
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	DLNSSTDSAA SASAASDISV ESKEATVCTV TLEKTSAGLG FSLEGGKGS L HGDKPLTINR IFKGTEQGEM VQPGDEILQL AGTAVQGLTR FEAWNVIKAL PDGPVTIVIR RTSLQCKQTT ASADS (125)
Purity	98%
Source	Yeast
MW	13.0 kDa
Entrez	16170
Expiration Date	6 months from date of receipt.