

Product Datasheet

Mouse IL-16 Recombinant Protein (orb1819790)

| | |
|-------------------------|--|
| Description | The Mouse IL-16 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The... |
| Target | IL-16 |
| Form/Appearance | Lyophilized |
| Storage | -20°C |
| Note | For research use only |
| Protein Sequence | DLNSSTDSAA SASAASDISV ESKEATVCTV TLEKTSAGLG FSLEGGKGS L HGDKPLTINR IFKGTEQGEM VQPGDEILQL AGTAVQGLTR FEAWNVIKAL PDGPVTIVIR RTSLQCKQTT ASADS (125) |
| Purity | 98% |
| Source | Yeast |
| MW | 13.0 kDa |
| Entrez | 16170 |
| Expiration Date | 6 months from date of receipt. |