

Product Datasheet

Zebrafish TGF-alpha Recombinant Protein (orb1819793)

| | |
|-------------------------|---|
| Description | The Zebrafish TGF-alpha yeast-derived recombinant protein is not tagged and is naturally endotoxin-f... |
| Target | TGF alpha |
| Form/Appearance | Lyophilized |
| Storage | -20°C |
| Note | For research use only |
| Protein Sequence | SHFDDCPDSH SHFCFHGTCR FLILEETPAC VCHPGFVGMR CEHADLLA (48) |
| Purity | 98% |
| Source | Yeast |
| MW | 5.4 kDa |
| Entrez | 100034402 |
| Expiration Date | 6 months from date of receipt. |