

## Product Datasheet

# Atlantic Salmon TGF-alpha Recombinant Protein (orb1819794)

<b>Description</b>	The Atlantic Salmon TGF-alpha yeast-derived recombinant protein is not tagged and is naturally endot...
<b>Target</b>	TGF alpha
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	VRSHFDDCPD SHSHFCFHGT CRFLILEETP ACVCHQGFVG MRCEHADLLA (50)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	5.7 kDa
<b>Entrez</b>	<b>106583550</b>
<b>Expiration Date</b>	6 months from date of receipt.