

Product Datasheet

Atlantic Salmon TGF-alpha Recombinant Protein (orb1819796)

Description	The Atlantic Salmon TGF-alpha yeast-derived recombinant protein is not tagged and is naturally endot...
Target	TGF alpha
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	VRSHFDDCPD SHSHFCFHGT CRFLILEETP ACVCHQGFVG MRCEHADLLA (50)
Purity	98%
Source	Yeast
MW	5.7 kDa
Entrez	106583550
Expiration Date	6 months from date of receipt.