

Product Datasheet

Human IL-33 Recombinant Protein (orb1819808)

| | |
|-------------------------|--|
| Description | The Human IL-33 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The... |
| Target | IL-33 |
| Form/Appearance | Lyophilized |
| Storage | -20°C |
| Note | For research use only |
| Protein Sequence | SITGISPITE YLASLSTYND QSITFALEDE SYEIIYVEDLK KDEKKDKVLL SYYESQHPSN ESGDGVDGKM LMVTLSPKTD FWLHANNKEH SVELHKCEKP LPDQAFFVLH NMHSNCVSFE CKTDPGVFIG VKDNHLALIK VDSSENLCCTE NILFKLSET (159) |
| Purity | 98% |
| Source | Yeast |
| MW | 18.0 kDa |
| Entrez | 90865 |
| Expiration Date | 6 months from date of receipt. |