

Product Datasheet

Human IL-33 Recombinant Protein (orb1819809)

Description	The Human IL-33 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	IL-33
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SITGISPITE YLASLSTYND QSITFALEDE SYEIIYVEDLK KDEKKDKVLL SYYESQHPSN ESGDGVDGKM LMVTLSPKTD FWLHANNKEH SVELHKCEKP LPDQAFFVLH NMHSNCVSFE CKTDPGVFIG VKDNHLALIK VDSSENLCCTE NILFKLSET (159)
Purity	98%
Source	Yeast
MW	18.0 kDa
Entrez	90865
Expiration Date	6 months from date of receipt.