

## Product Datasheet

### Human IL-33 Recombinant Protein (orb1819810)

<b>Description</b>	The Human IL-33 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
<b>Target</b>	IL-33
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	SITGISPITE YLASLSTYND QSITFALEDE SYEIIYVEDLK KDEKKDKVLL SYYESQHPSN ESGDGVDGKM LMVTLSPKTD FWLHANNKEH SVELHKCEKP LPDQAFFVLH NMHSNCVSFE CKTDPGVFIG VKDNHLALIK VDSSENLCTE NILFKLSET (159)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	18.0 kDa
<b>Entrez</b>	<b>90865</b>
<b>Expiration Date</b>	6 months from date of receipt.