

## Product Datasheet

# Rat CD40 Ligand (CD40LG) Recombinant Protein (orb1819814)

<b>Description</b>	The Rat CD40 Ligand yeast-derived recombinant protein is not tagged and is naturally endotoxin-free....
<b>Target</b>	CD40LG
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	MQKGDEDPQI AAHVVSEANS NAASVLQWAK KGYTMTKSNL VLENGRQLT VKREGLYYVY TQVTFCSNRE PLSQRPFIVS LWLKPSSGSE RILLRAANTH SSSKLCEQQS IHLGGVFELQ AGASVFNVT EASQVIHGIG FSSFGLLKL (149)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	16.4 kDa
<b>Entrez</b>	<a href="#">84349</a>
<b>Expiration Date</b>	6 months from date of receipt.