

Product Datasheet

Rabbit TGF-alpha Recombinant Protein (orb1819818)

Description	The Rabbit TGF-alpha yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
Target	TGF alpha
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	VVSHFNQCPD SHTQFCFHGT CRFLVQEDKP ACVCHSGYVG ARCEHADLLA (50)
Purity	98%
Source	Yeast
MW	5.6 kDa
Expiration Date	6 months from date of receipt.