

Product Datasheet

Human IL-31 Recombinant Protein (orb1819829)

Description	The Human IL-31 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	IL-31
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SHTLPVRLLR PSDDVQKIVE ELQSLSKMLL KDVEEEKGVL VSQNYTLPLCL SPDAQPPNNI HSPAIRAYLK TIRQLDNKSV IDEIIEHLDK LIFQDAPETN ISVPTDTHC KRFILTISQQ FSECMDLALK SLTSGAQQAT T (141)
Purity	98%
Source	Yeast
MW	15.8 kDa
Entrez	386653
Expiration Date	6 months from date of receipt.