

Product Datasheet

Ovine CD40 Ligand (CD40LG) Recombinant Protein (orb1819833)

Description	The Ovine CD40 Ligand yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
Target	CD40LG
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	MHKGDQEPQI AAHVISEASS KTTSVLQWAP KGYTTLSSNL VTLENGKQLA VKRQGFYYIY TQVTFCSNRE ALNQAPFIAS LCLKSPSGSE RILLRAANTH SSSKPCGQQS IHLGGVFELQ QGASVFNVT DPSQVSHGTG FTSFGLLKL (149)
Purity	98%
Source	Yeast
MW	16.1 kDa
Entrez	767628
Expiration Date	6 months from date of receipt.