

Product Datasheet

Ferret Midkine Recombinant Protein (orb1819839)

Description	The Ferret Midkine yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. ...
Target	Midkine
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	KKKDKVKKGG PGTECAEWAW GPCTPSSKDC GVGFRGTCG AQTRTRCRV PCNWKKEFGA DCKYKFESWG ACDGGTGTKA RQGLKARY NAQCQETIRV TKPCTPKTKA KAKAKKGGK D (121)
Purity	98%
Source	Yeast
MW	13.2 kDa
Entrez	101682417
Expiration Date	6 months from date of receipt.