

Product Datasheet

Canine, Equine, Feline, Ferret, and Swine TGF-alpha Recombinant Protein (orb1819843)

| | |
|-------------------------|---|
| Description | The Canine, Equine, Feline, Ferret, and Swine TGF-alpha yeast-derived recombinant protein is not tag... |
| Target | TGF alpha |
| Form/Appearance | Lyophilized |
| Storage | -20°C |
| Note | For research use only |
| Protein Sequence | VVSHFNDCPD SHSQCFHGT CRFLVQEDKP ACVCHSGYVG ARCEHADLLA (50) |
| Purity | 98% |
| Source | Yeast |
| MW | 5.5 kDa |
| Expiration Date | 6 months from date of receipt. |