

Product Datasheet

Rat IL-1 alpha Recombinant Protein (orb1819847)

| | |
|-------------------------|--|
| Description | The Rat IL-1 alpha yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. ... |
| Target | IL-1 alpha |
| Form/Appearance | Lyophilized |
| Storage | -20°C |
| Note | For research use only |
| Protein Sequence | SAPHSFQNNL RYKLIRIVKQ EFIMNDSLNQ NIYVDMDRIH LKAASLNDLQ LEVKFDMYAY SSGGDDSKYP VTLKVSNTQL FVSAQGEDKP VLLKEIPETP KLITGSETDL IFFWEKINSK NYFTSAAFPE LLIATKEQSQ VHLARGLPSM IDFQIS (156) |
| Purity | 98% |
| Source | Yeast |
| MW | 17.8 kDa |
| Entrez | 24493 |
| Expiration Date | 6 months from date of receipt. |