

Product Datasheet

Feline IFN alpha Recombinant Protein (orb1819852)

Description	The Feline IFN alpha yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
Target	IFN alpha
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	CDLPQTHGLL NRRALTLGQ MRRLPASSCQ KDRNDFAFPQ DVFGGDQSHK AQALSVVHVT NQKIFHFFCT EASSSAAWNT TLLEEFCTGL DRQLTRLEAC VLQEVEEGEA PLTNEDIHPE DSILRNYFQR LSLYLQEKKY SPCAWIVRA EIMRSLYSS TALQKRLRSE K (171)
Purity	98%
Source	Yeast
MW	19.7 kDa
Entrez	493653
Expiration Date	6 months from date of receipt.