

Product Datasheet

Human IFN alpha Recombinant Protein (orb1819856)

Description	The Human IFN alpha yeast-derived recombinant protein is not tagged and is naturally endotoxin-free....
Target	IFN alpha
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	CDLPETHSLD NRRTLMLLAQ MSRISPSSCL MDRHDFGFPPQ EEFDGNQFQK APAISVLHEL IQQIFNLFTT KDSSAAWDED LLDKFCTELY QQLNDLEACV MQEERVGETP LMNADSILAV KKYFRRITLY LTEKKYSPCA WEVVRAEIMR SLSLSTNLQE RLRRKE (166)
Purity	98%
Source	Yeast
MW	19.4 kDa
Entrez	3447
Expiration Date	6 months from date of receipt.