

Product Datasheet

Dolphin CCL4 Recombinant Protein (orb1819859)

Description	The Dolphin CCL4 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	CCL4
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	LEKRAPMSSD PPTACCFYTYT LRKLPRHFVI DYYETSSLCS QPAVVFQTKK GRQVCANPSD PWVQEYMDDL ELN (73)
Purity	98%
Source	Yeast
MW	8.4 kDa
Entrez	101330472
Expiration Date	6 months from date of receipt.