

Product Datasheet

Human M-CSF Recombinant Protein (orb1819863)

Description	The Human M-CSF yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	M-CSF
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	EEVSEYCSHM IGSGHLQSLQ RLIDSQMETS CQITFEFVDQ EQLKDPVCYL KKAFLLVQDI MEDTMRFRDN TPNAIAIVQL QELSLRLKSC FTKDYEEHDK ACVRTFYETP LQLEKVKNV FNETKNLLDK DWNIFSKNCN NSFAECSSQG HERQSEGS (158)
Purity	98%
Source	Yeast
MW	18.4 kDa
Entrez	1435
Expiration Date	6 months from date of receipt.