

Product Datasheet

Dolphin CCL14 Recombinant Protein (orb1819866)

Description	The Dolphin CCL14 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. T...
Target	CCL14
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SKTESDSRGP YPTECCYSY LTRPIRKHWI SRYYETGGLC SKPGIVFITK NGRKLCANPS DDWVQDYIKG AGGDQSDPEV AEKDTAQNLIK RWG (93)
Purity	98%
Source	Yeast
MW	10.5 kDa
Entrez	101338296
Expiration Date	6 months from date of receipt.