

Product Datasheet

Ovine/Caprine Amphiregulin Recombinant Protein (orb1819874)

Description	The Ovine & Caprine Amphiregulin yeast-derived recombinant protein is not tagged and is naturally en...
Target	Amphiregulin
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SVRVEQVVKP KKNKTESEKT SDKPKRKKKG GKNGKNRRNR KRKNLCDTEF QNFCIHGKCT FLEQLETVSC QCYPEYFGER CGEKSMK (87)
Purity	98%
Source	Yeast
MW	10.2 kDa
Entrez	102180661 , 443427
Expiration Date	6 months from date of receipt.