

Product Datasheet

Swine Amphiregulin Recombinant Protein (orb1819877)

Description	The Swine Amphiregulin yeast-derived recombinant protein is not tagged and is naturally endotoxin-fr...
Target	Amphiregulin
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SVRVEQVVKP KRNRTESNT SDKPKRKKKG GKSGKNRRNR KKKNPCDAEF QNFCIHGDCK YIEHLEAVTC KCYQDYFGER CGEKSMK (87)
Purity	98%
Source	Yeast
MW	10.2 kDa
Entrez	397668
Expiration Date	6 months from date of receipt.