

Product Datasheet

Ferret TNF alpha Recombinant Protein (orb1819887)

Description	The Ferret TNF alpha yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
Target	TNF alpha
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	LQWLSRRANA LLANGVELTD NQLIVPSDGL YLIYSQVLFK GRGCSSTNVL LTHTISRFAV SYRTRKVNLLS AIKSPCQRET PEGTEAKPWY EPIHLGGVFQ LEKGDRLSAE INLPAYLDFA ESGQVYFGII AL (132)
Purity	98%
Source	Yeast
MW	14.6 kDa
Entrez	101687148
Expiration Date	6 months from date of receipt.