

## Product Datasheet

### Feline CCL3 Recombinant Protein (orb1819888)

<b>Description</b>	The Feline CCL3 Amphiregulin yeast-derived recombinant protein is not tagged and is naturally endoto...
<b>Target</b>	CCL3
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	APFGADTPTA CCFSYVSKQI PRKFIADYFE TSSQCSKPGV IFQTRRGRQV CADPSEAWVQ EYVTDLELSA (70)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	7.8 kDa
<b>Entrez</b>	<a href="#">100302540</a>
<b>Expiration Date</b>	6 months from date of receipt.