

Product Datasheet

Cynomolgus Monkey Amphiregulin Recombinant Protein (orb1819891)

Description	The Cynomolgus Monkey Amphiregulin yeast-derived recombinant protein is not tagged and is naturally ...
Target	Amphiregulin
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SVRVEQVVKP PKNKTESENT SDKPKRKKKG GKNGKNRRNR KKKNPCNAEF QNFCIHGECK YIEHLEAVTC KCQQEYFGER CGEKSMK (87)
Purity	98%
Source	Yeast
MW	10.1 kDa
Entrez	102118341
Expiration Date	6 months from date of receipt.