

Product Datasheet

Rabbit Amphiregulin Recombinant Protein (orb1819892)

Description	The Rabbit Amphiregulin yeast-derived recombinant protein is not tagged and is naturally endotoxin-f...
Target	Amphiregulin
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SVRVEQVIKP KKNKAESEKT SDKPKRKKKG GKNGKNRKNK KKKNPCDAEF QNFCIHGECK YIEHLDAVTC KCHHDYFGER CGEKSMK (87)
Purity	98%
Source	Yeast
MW	10.1 kDa
Entrez	100354550
Expiration Date	6 months from date of receipt.