

Product Datasheet

Cynomolgus Monkey CCL3 Recombinant Protein (orb1819905)

Description	The Cynomolgus Monkey CCL3 yeast-derived recombinant protein is not tagged and is naturally endotoxi...
Target	CCL3
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	NRISATFAAD TPTSCCFSYI SRQIPQNFIA DYFETNSQCS KPGVIFLTR GRQVCADPSK EWWQKYVSDL ELSA (74)
Purity	98%
Source	Yeast
MW	8.3 kDa
Entrez	102136134
Expiration Date	6 months from date of receipt.