

Product Datasheet

Feline Amphiregulin Recombinant Protein (orb1819909)

Description	The Feline Amphiregulin yeast-derived recombinant protein is not tagged and is naturally endotoxin-f...
Target	Amphiregulin
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SVRVEQVVKP MKNKTESEKT SDKPKRKKKG GKNGKNRKNR KKKNPCDAEF QNFCIHGECK YIENLGAVTC KCHQDYFGER CGEKSMK (87)
Purity	98%
Source	Yeast
MW	10.0 kDa
Entrez	101081041
Expiration Date	6 months from date of receipt.